## Bebe-Ball



Bebe-Ball Manual Suction Device		
Order #	Product Description	Bx Qty
Bebe-Ball with Bebeonkers		
16-4	Bebe-Ball with Bebeonkers 16-2 (Std)	50
16-5	Bebe-Ball with UP Bebeonkers 16-1 (bulb)	50
16-6	Bebe-Ball with Webeonkers 16-0 (UP)	50
Bebe-Ball with Bebeonkers and Keeper Kover		
16-4KK	Bebe-Ball with Bebeonkers 16-2 (Std) with Keeper Kover	50
16-5KK	Bebe-Ball with UP Bebeonkers 16-1 (bulb) with Keeper Kover	50
16-6KK	Bebe-Ball with Webeonkers 16-0 (UP) with Keeper Kover	50



## Bebe-Ball


- Bulb suction adaptor for Bebeonkers 16-0, 16-1 and 16-2
- Helps prevent damage that occurs if one-hole device accidently touches the sensitive mucosa of a newborn or preemie
- Disposal after 12 to 24 hours or per hospital policy to prevent infection is recommended
- Use on transports or when power is unavailable and STAT suctioning is required
- Standard bulb syringe has one hole and is inappropriate because it can grab and tear tender tissue too



Red ball is our trial color, actual production will be blue translucent.

## Bebe Keeper Kover

Can be used up to 12 hours on the same patient, or per hospital policy.




**Bx Qty** 

300



Bebe Keeper Kover is an effective cover to keep the Bebeonkers clean.

Can be used on Bebeonkers I and Bebeonkers II with suction on or off, protecting the device while not in use.



Can also be used with the Bebe-Ball.

Note new blue translucent color